

MVE CryoCart



Features include

- Temperature monitor with type T Thermocouple
- Optional Temperature monitor and LN₂ transfer hose connection with plumbing system (manual lever)

MVE CryoCarts are designed for the loading of biological samples into canes, boxes, racks or frames. When used as a portable workbench, the unit will provide a safe and controlled environment for your samples for up to eight hours with the lid open. MVE CryoCart is ideal for transporting large quantities of samples from one tank to another within the same facility.



MVE CryoCart	
Performance	
LN ₂ Capacity w/o Inventory To Top of Platform L. est.	39.5
Unit Dimensions	
Overall Length in. (mm)	54.46 (1408)
Overall Width in. (mm)	20.16 (512)
Overall Height in. (mm)	39.5 (1003)
Lift Overall Height in. (mm)	37.2 (944)
Inside Length in. (mm)	41.0 (1041)
Inside Width in. (mm)	13.2 (335)
Inside Depth in. (mm)	15.4 (391)
Platform Height in. (mm)	4.5 (114)
Footprint in. (mm)	17.6x51.7 (447x1313)
Weight Empty lb. (kg) est.	195 (88)
Weight w/LN ₂ Filled to Top of Platform lb. (kg) est.	266 (121)
Number of Racks (15-2)	2

**90 Days Parts Warranty • One Year Vacuum Warranty
Not designed for long-term storage.**

www.mvebio.com
844-MVE-CRYO